

Products for Dentistry

OSUNG Catalogue 2017/2018

**Instrument
Management**


Instrument Cassettes

1. No lid
2. Fix with pressure bar
3. Easy to clean and set instruments
4. Silicone on bottom prevent damage of products

Economical & essential product for user of wash machine & ultrasonic cleaner

EFS8

- Exterior dimensions : 140 x 187 x 25H(mm)
- 8 instruments capacity of periodontal & diagnostic (Max. length 184mm)




Essential for protection of blades. The best cost-effectiveness product

EFS8L

- Exterior dimensions : 222 x 140 x 25H(mm)
- 8 instruments capacity (Max. length 221mm)

Charactor

Designed for long-sized instruments like periosteal elevators, sinus lifts and so on.




Simple type & high economics

EFS16

- Exterior dimensions : 281 x 187 x 25H(mm)
- 16 instruments capacity of periodontal & diagnostic (Max. length 184mm)

Charactor

Possible to make implant surgical kit with osteotome if needed



45JD

50JD

65 Jd

Instrument Cassettes



1. Silicone on bottom prevent damage of products.
2. Spring lock system is easy to open & close.


For students of dental hygienic

EFCCN1

- Exterior dimensions : 139 x 185 x 21.7H(mm)
- 10 periodontal instruments capacity

Charactor

It comes in with a lid and is compact, easy to keep, and portable. This product is made of stainless steel with much consideration given to structure and shape which can sustain big loads and impacts for students to make convenient use of it.



For composition of surgical operation kit enhanced durability

EFCCCL1

- Exterior dimensions : 305 x 210 x 34H(mm)

Charactor

A cassette for periodontal and implant surgery. The lock employs a spring lock system, which offers easy opening and closing without mechanical failures.




For composition of surgical operation kit

Pls check whether it is suitable to your autoclave as it is large size's cassette.

EFCCCL2

- Exterior dimensions : 375 x 240 x 34H(mm)



65JD

110JD


125JD

Instrument Cassettes

This is a cassette designed to store such large tools as forceps and elevators

EFCCL1-F

- Exterior dimensions : 305 x 210 x 34H(mm)



Charactor

1. For surgical operation
2. Silicone pad insert for protection of product
3. Spring lock type locking system


This cassette is designed for the storage of elevators and luxators

EFCCL1-L

- Exterior dimensions : 305 x 210 x 34H(mm)

Charactor

1. For surgical operation
2. Silicone pad insert for protection of product
3. Spring lock type locking system



Cassette specifically designed to optimize cleaning, sterilization and instrument protection

NEW PRODUCT

EFCCL15

- Exterior dimensions : 205 x 271 x 34H(mm)
- 15 instruments capacity of periodontal & diagnostic

Charactor

1. Easy to clean instruments
2. Silicone pad insert for protection of product



110JD

110JD

110JD

Surgi-Drill Stand

Surgi-Drill Stand


Surgi-Drill Stand is a perfect solution for managing surgical drills and burs. The stand has 16 multi-silicone-holders and the multi-holder can hold any kinds of drill & bur.

DSTA16

• Size 98 x 40 x 45H(mm)




The ong-sized drill can be contained by overturning the stand. Long drills can be contained like the picture, as the some part of the bottom is windowed for long drills. In this case, the upper side of the cover becomes the bottom.



35JD

1.



Multi-holder : The centered large hole is for a drill or a low speed bur, and the winged small holes are for high speed burs.

One holder can hold two FG burs.

Inserted implant tissue punch

2.



The cover can be used as a mini-tray when it is opened.

Surgical Drape · Wrapping Clothes · Instrument Pouches

50JD

Surgical Drape



- Convenient drape for surgery. A magnet is inserted inside of drape to stick hinge-typed instruments such as scissors, needle holder and so on.



WDMA

- Size : 900 x 900(mm)
- Hole diameter : 90(mm)

* When sterilization, pls take out magnet and put magnet again after autoclave.

Wrapping Cloth



- No discolor during washing. Made in Korea

WR5050

- Wrapping cloth for sterilizing dental instruments
- Size : 500 x 500(mm)

13JD

WR7575

- Wrapping cloth for sterilizing dental instruments
- Size : 750 x 750(mm)

16JD



Videp
Clip

Durable fabric and luxurious color! Best choice for those who prefer a cloth pouch for instrument storage and sterilization



- Used to store and sterilize dental instruments. Also it can be used as instrument cassette during operation. It is recommended to the dentists who prefer soft material to hard material like metal cassette. It is made of linen two layered with 14 embroidered.

WPA

- Size : 470 x 400 (mm)

18 JD

WPB

- Size : 550 x 400 (mm)

Endo Ruler· Endo File Boxes

Endo Bath

It holds up to 44 files and reamers for autoclave sterilization.

50JD



EAA1

Endo Bath, Silver
· Size 78 x 59(H)mm

Endo Can

It is a small sized box for each patient for preventing infection.
It holds up to 20 files and reamers for autoclave sterilization.

45JD



EAB1

Endo Can, Silver
· Size 46 x 60(H)mm


Scaler Tip Management

18 JD

Scaler Tip Stand

USSTA5

- Size : 50 x 40 x 38H(mm)
- Used for sterilization & storage of ultrasonic scaler tips. 5 tips capacity.



CHARACTER

1. Superior durability by stainless steel
2. Easy to move and sterilize by small size and use it on chair-side

Scaler Tip Torque Wrench

USETW

- EMS / SATELEC



18 JD

CHARACTER

1. Use as tip stand during sterilization
2. It is safe from infection as tip is not touch to hand when scaler tip connects to handpiece.
3. Handpiece is protected by connecting tip with regular power.

Sharpening Stone

NEW PRODUCT

SST-C3

- Ceramic Sharpening Stone #3C (Medium Grif)
- Brown
- 80 x 33 x 6.3H (mm)



50JD

Bur Block

EBB1

Bur Block, Silver

- Size 56 x 26 x 29(H)mm
- It can be shorten treatment time to prepare kits for each treatment as FG burs & RA burs can be kept together.
- [Autoclavable](#)



18 JD


Bur Management

Free Bur Stand

132°C Autoclavable

EBFS

- Size 300 x 60mm
- Color : ■ ■ ■ Blue, Red, Yellow



CHRACTOR

1. Capacity of 40 holes for FG & RA burs and possible to cut desired size.
2. Easy to cut with medical silicone and no loss of bur by elasticity of silicone.


Practice



Put various burs



Cutting bur sand



Inserted bur



No drop even the stand is overturned

Bur Block

132°C Autoclavable

NEW PRODUCT

EBSTAFR87

- Size 61 x 15 x 29H(mm)



Max height of bur 25.5mm



RA
FG

FG Bur 8ea & RA Bur 7ea

Merits of EBSTAFR87

- Make better use of space as compact size
- Capacity of FG Bur 8ea & RA Bur 7ea at a same time.
- No loss during keeping or moving by hinge type's cover
- Use as personal kit, it helps to prevent infection.
- Locking Standing

15JD

15JD

Cotton Pellet Devices

Cotton Pellet Dispenser

RGCPD

- Size 45 x 60 x 54H(mm)



35JD

Cotton Pellet Push Device

- One action convenient to use

RGCP1

- Stainless Steel
- Autoclavable
- Size 45 x 45 x 48H(mm)



15 JD

One touch! Cotton Pellet Push Device pushes up cotton pellet automatically using cotton pellet dispenser's weight so it is easy to pick cotton pellet up by pincette.

