




WOMEN'S HEALTH


INTRODUCTION	2
OVUM PICKUP NEEDLES	
ACE Single Lumen Ovum Pickup Needle.....	3
ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line.....	4
ACE-M Manual Single Lumen Ovum Pickup Needle.....	5
OVUM PICK Double Lumen Ovum Pickup Needle.....	6
EMBRYO TRANSFER CATHETERS	
EMBRYO TRANS Embryo Transfer Catheter.....	7
ECHO TRANS Embryo Transfer Catheter.....	8
TWINKLE EMBRYO TRANS Embryo Transfer Catheter.....	9
BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip.....	10
BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip.....	11
PRO EMBRYO TRANS Embryo Transfer Catheter.....	12
PRO ECHO TRANS Embryo Transfer Catheter.....	13
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter.....	14
INTRA UTERINE INSEMINATION CATHETERS	
INSEMINA Intra Uterine Insemination Catheter.....	15
IUI CURVED Intra Uterine Insemination Catheter Curved	16
TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet.....	17
FLEXI Intra Uterine Insemination Catheter.....	18
STUMPY Intra Uterine Insemination Cannula.....	19
FALLOPIAN TUBE RECANALISATION CATHETERS	
SALPINX CATH Fallopian Tube Recanalisation Catheter.....	20
HYSTEOSALPINGOGRAPHY CATHETER	
HYSTERO CATH Hystero Salpingography Catheter.....	21
ENDOMETRIAL BIOPSY SAMPLING SYSTEM	
THE SAMPLER Endometrial Biopsy Catheter.....	22
PRO-VAC Endometrial Biopsy Sampler.....	23
VITRIFICATION STRAWS	
SURESTRIP Vitrification Straws (Open System Vitrification).....	24
SURELOCK Vitrification Straws.....	25
SUREHOOK Vitrification Straws.....	26
OTHER PRODUCTS	
ASPI Aspiration Needle.....	27

ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line




Allwin ovum pickup needles are used for laparoscopic or ultrasound guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line
- Alternate aspiration & flushing made possible through single channel
- Multiple operation functionality by using stopcock valve




Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	16	10
A-OVM-SL17/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	17	10
A-OVM-SL18/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	18	10

*All products may not be available in all countries. Please contact our sales department for details.




Highly echogenic tip

ECHO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Soft inner catheter with smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Visible under ultrasound
- Stylet for outer sheath available separately



Soft inner catheter with smooth tip

Marking on inner catheter and outer sheath


Stylet for Embryo Trans / Echo Trans Catheter / Twinkle Embryo Trans

Order No.	Description	Catheter Length	Qty. / Box
ET-ETC23	ECHO TRANS Embryo Transfer Catheter	23	10
ET-ETC18	ECHO TRANS Embryo Transfer Catheter	18	10
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10

*All products may not be available in all countries. Please contact our sales department for details.



PRO ECHO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

Allwin embryo transfer catheters (pro embryo trans, pro echo trans, and twinkle pro echo trans) also include stylets or soft obturators intended to assist uterine access of the guide catheter during an embryo transfer procedure.

- Soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible transfer catheter
- Visible under ultrasound
- Orientation arrow on hub



Soft obturator for Pro Embryo Trans / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans




Stylet for Pro Embryo Transfer / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans

Order No.	Description	Qty. / Box
PET-ETC	PRO ECHO TRANS Embryo Transfer Catheter	10
PET-ETC-SB	PRO ECHO TRANS Embryo Transfer Catheter with Soft Obturator	10
PET-ETC-STY	PRO ECHO TRANS Embryo Transfer Catheter with Stylet	10
SB-PET-ETC	Soft obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Transfer / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.




TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible transfer catheter
- Silicone depth positioner for precise placement
- Orientation arrow on hub



Soft obturator for Pro Embryo Trans / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans Catheter



Stylet for Pro Embryo Trans / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans Catheter



Order No.	Description	Qty. / Box
TPEB-ETC	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter	10
TPEB-ETC-SB	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
TPEB-ETC-STY	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Transfer / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.



allwin[®]

Medical Devices



US Sales

✉ sales.usa@allwinmedical.com

☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com

☎ +91 9987084563

Worldwide sales (Excluding India/US)

✉ sales@allwinmedical.com

☎ +91 95949-81581