



International Quality IUI Media From Shivani Scientific

Safe, Potential,
Efficient and Convenient
approach for
Sperm Preparation

Dedicated cold chain & 24 hours
temperature controlled supply

Please register for IUI training

- Media is manufactured at environmentally controlled, ISO9001 & ISO13485 certified facility in EUROPE.
- Repacked and labelled at our World class State-of-art, ART Medical Devices Facility at Mira, Mumbai-104, India.
- The Manufacturing facility comply CDSCO (FDA) and international standard, one of its kind in INDIA.





SPERM SCREEN™ KIT

Ready-to-use 3 in 1 kits for gradient Sperm Separation

- Sperm Screen is an isotonic gradient for semen preparation with silane-coated colloidal silica particles suspended in HEPES Buffered Isotonic Salt Solution.
- With distinct colour coding of labels for identification. ●●●
- SPERM SCREEN Kit has 10 kits.



SPERM PRIME™

- Sperm Prime is sterile balanced salt solution, HEPES buffered semen wash medium of Human Serum Albumin (hSA)
- Optimised for 'Sperm Wash' and medium for preparing sperm.
- Ensure effortless pipetting, Smooth aspiration and Warming at 37 C.